

The Inman Orthodontic Family of 3-3'S

The UK Style 3-3

The UK Style 3-3 is our most popular design.

The jig holds the flat braided soft wire in place so the clinician can add composite to each tooth, from 3 to 3. The jig is then cut away and the composite is finished and polished smooth.



The Direct Bond 3-3

.030 stainless round wire with bonding pads on the cuspids. This appliance is usually bonded to the cuspids (as shown). Fabricated with or without a putty transfer jig.



If fabricated without the bonding pads, helices or a right angle bend can be used in place of the bonding pads. .028 TMA wire may also be used.

Your trusted source for orthodontic appliances.



inmanortho.com

1800.289.0118

Inman Orthodontic Laboratory Inc
3953 NW 126th Ave
Coral Springs, FL 33065

The Inman 3-3's

The Rectangular Stiff 3-3

Fabricated with a stiff rectangular wire.

FEATURES:

With or without a putty jig.

Composite is added to all the linguals, cuspid to cuspid.

Some clinicians use floss pulled through the interproximal as a jig.



The Hygienic 3-3

Light round wire contoured to allow easy flossing of the interproximals.

Fabricated with the putty jig. The cuspids are bonded, then the jig is removed.



The Printed 3-3

This is a one piece **CAD/CAM Cobalt** printed appliance.

A bonding pad is placed on the lingual of all the anterior teeth.

A putty jig can be used but due to the intimate fit it can easily be placed directly.



Your trusted source for orthodontic appliances.



inmanortho.com

1800.289.0118

Inman Orthodontic Laboratory Inc
3953 NW 126th Ave
Coral Springs, FL 33065