



COSMO IPC

INTRODUCING THE COSMO IPC

A totally new approach to removable expansion appliances.

The gentle, steady force slowly and efficiently expands the arch without creating discomfort to the patient. The Cosmo IPC is ideal for patients in the early/mixed dentition, revolutionizing maxillary skeletal and dental development in a predictable noninvasive manner.

In development for over 5 years and proven in clinical trials, the Cosmo IPC uses open coil nickel titanium spring to generate orthodontic force.

The appliance is adjusted by the clinician eliminating patient or parent obligations.

Case studies are available @ inmanortho.com

The CosmoIPC is fabricated with the award winning Inman Power Component.

The IPC's sliding lock compresses open coil spring to generate gentle orthodontic force.

Designed with the guidance of Dr. Cosmo Haralambidis.



inmanortho.com
1800.289.0118

Your trusted source for orthodontic appliances