


# INTRODUCING THE **IDS** Inman Digital Splint

Our new ID Splint is fabricated with a completely digital process. This eliminates the shrinkage and warping issues that cause splints to rock and not fit properly. The ID Splint is digitally designed and digitally printed to your prescription and the fit is absolutely perfect!



## Inman Digital Splint Precise. Durable. Effective.

- **ID Splints** are computer designed and digitally manufactured.
- **INTRAORAL SCAN** compatible. Upload your scans to our cloud based Sharefile.
- **NO SCANNER**, no problem. We will digitize your models and use your construction bite.

## All popular upper and lower splint designs supported.

We mount the digital models on a virtual articulator and design the ID Splint to your RX. Splints are returned on 3D - printed models.

Lost splints can be reordered by phone and remade from archived files at a discount.



Call Inman Orthodontic Laboratory to learn more.



Experience the difference...  
your patients will appreciate the results!

**1.800.289.0118**

[www.inmanortho.com](http://www.inmanortho.com)  
[www.inmandigner.com](http://www.inmandigner.com)